



tasarimmed[®]



EXTERNAL CIRCULAR FIXATION SYSTEM (ADULT)

EXTERNAL CIRCULAR FIXATION SYSTEM (ADULT)



HALF RINGS CARBON FIBRE

YH-4050 Ø 50 MM	YH-4150 Ø 150 MM
YH-4060 Ø 60 MM	YH-4160 Ø 160 MM
YH-4080 Ø 80 MM	YH-4180 Ø 180 MM
YH-4100 Ø 100 MM	YH-4200 Ø 200 MM
YH-4120 Ø 120 MM	YH-4220 Ø 220 MM
YH-4130 Ø 130 MM	YH-4240 Ø 240 MM
YH-4140 Ø 140 MM	



5/8" C - RINGS CARBON FIBRE

CH-4025 Ø 25 MM	CH-4130 Ø 130 MM
CH-4030 Ø 30 MM	CH-4140 Ø 140 MM
CH-4035 Ø 35 MM	CH-4150 Ø 150 MM
CH-4045 Ø 45 MM	CH-4160 Ø 160 MM
CH-4050 Ø 50 MM	CH-4180 Ø 180 MM
CH-4060 Ø 60 MM	CH-4200 Ø 200 MM
CH-4080 Ø 80 MM	CH-4220 Ø 220 MM
CH-4100 Ø 100 MM	CH-4240 Ø 240 MM
CH-4120 Ø 120 MM	



RING WITH CURVED EXTREMITIES (FOOT RING) CARBON FIBRE

OH-4080 Ø 80 MM
OH-4100 Ø 100 MM
OH-4120 Ø 120 MM
OH-4130 Ø 130 MM
OH-4140 Ø 140 MM
OH-4150 Ø 150 MM



FA-4000 FEMURAL ARCH CARBON FIBRE



HALF RING SPECIAL AL. ALLOY

YH-3025 Ø 25 MM	YH-3120 Ø 120 MM
YH-3030 Ø 30 MM	YH-3130 Ø 130 MM
YH-3035 Ø 35 MM	YH-3140 Ø 140 MM
YH-3040 Ø 40 MM	YH-3150 Ø 150 MM
YH-3045 Ø 45 MM	YH-3160 Ø 160 MM
YH-3050 Ø 50 MM	YH-3180 Ø 180 MM
YH-3060 Ø 60 MM	YH-3200 Ø 200 MM
YH-3080 Ø 80 MM	YH-3220 Ø 220 MM
YH-3100 Ø 100 MM	YH-3240 Ø 240 MM



5/8" C - RINGS SPECIAL AL. ALLOY

CH-3025 Ø 25 MM	CH-3130 Ø 130 MM
CH-3030 Ø 30 MM	CH-3140 Ø 140 MM
CH-3035 Ø 35 MM	CH-3150 Ø 150 MM
CH-3040 Ø 40 MM	CH-3160 Ø 160 MM
CH-3045 Ø 45 MM	CH-3180 Ø 180 MM
CH-3050 Ø 50 MM	CH-3200 Ø 200 MM
CH-3080 Ø 80 MM	CH-3220 Ø 220 MM
CH-3100 Ø 100 MM	CH-3240 Ø 240 MM
CH-3120 Ø 120 MM	



RING WITH CURVED EXTREMITIES (FOOT RING)

OH-3080 Ø 80 MM
OH-3100 Ø 100 MM
OH-3120 Ø 120 MM
OH-3130 Ø 130 MM
OH-3140 Ø 140 MM
OH-3150 Ø 150 MM



RING WITH CURVED EXTREMITIES (OMEGA RING)

OH-4000S SMALL
OH-4000M MEDIUM
OH-4000L LARGE

EXTERNAL CIRCULAR FIXATION SYSTEM (ADULT)



FEMORAL ARC

FA-3045	45°200 MM
FA-3090	90°200 MM
FA-3590	90°250 MM
FA-3120	120°200 MM
FA-3520	120°250 MM




OH-3001 EXTRA LARGE PELVIC C RING



CANNULATED THREADED RODS

DR-1060	L 60 MM
DR-1080	L 80 MM
DR-1100	L100 MM
DR-1120	L120 MM
DR-1150	L150 MM
DR-1200	L200 MM




THREADED RODS

DR-0060	L 60 MM
DR-0080	L 80 MM
DR-0100	L100 MM
DR-0120	L120 MM
DR-0130	L130 MM
DR-0150	L150 MM
DR-0200	L200 MM
DR-0300	L300 MM
DR-0400	L400 MM




SLOTTED THREADED RODS

OR-0040	L 40 MM
OR-0050	L 50 MM
OR-0060	L 60 MM
OR-0080	L 80 MM
OR-0100	L100 MM
OR-0150	L150 MM



GRADUATED TELESCOPIC RODS

TR-1100	L 100 MM
TR-1150	L 150 MM
TR-1200	L200 MM
TR-1250	L 250 MM



BP-0004 OBLIQUE SUPPORT CONNECTOR



BP-5000 EXTENSION PIECE



BP-5001 EXTENSION PIECE CONNECTION PLATE



KS-0400 QUADRAGONAL NUT



TF - 9000 HYBRID SLIDER CLAMP




BP-0050 SPECIAL QUADRAGONAL NUT



CONNECTION PLATES

BP-0038	2 HLS.	38 MM	BP-0104	8 HLS.	104 MM
BP-0049	3 HLS.	49 MM	BP-0115	9 HLS.	115 MM
BP-0060	4 HLS.	60 MM	BP-0126	10 HLS.	126 MM
BP-0071	5 HLS.	71 MM	BP-0159	8 HLS.	159 MM
BP-0082	6 HLS.	82 MM	BP-0236	12 HLS.	236 MM
BP-0093	7 HLS.	93 MM	BP-0336	17 HLS.	336 MM



CONNECTION PLATES WITH THREADED END

BP-0131	5 HLS.	131 MM
BP-0171	7 HLS.	171 MM
BP-0211	9 HLS.	211 MM
BP-0251	11 HLS.	251 MM



TWISTED PLATES

BP-0046	2 HLS.	46 MM
BP-0067	3 HLS.	67 MM
BP-0088	4 HLS.	88 MM



KE-1120 TWISTED HINGE


EXTERNAL CIRCULAR FIXATION SYSTEM (ADULT)



KE-0111
FEMALE HALF HINGE




KE-0100
MALE HALF HINGE



FEMALE SUPPORTS

KE-1001 1 HOLE
KE-1002 2 HOLES
KE-1003 3 HOLES
KE-1004 4 HOLES




MALE SUPPORTS

KE-0001 1 HOLE
KE-0002 2 HOLES
KE-0003 3 HOLES
KE-0004 4 HOLES



KY-0900
TRANSLATION
ROTATION BLOCK



KY-0001
MULTIPLE PIN
FIXATION BOLT



KY-0003
UNIVERSAL PIN
FIXATION BOLT



KY-0002
SINGLE PIN
FIXATION BOLT



PS-0016
WASHER
Ø 16 MM



KY-0005
ROUND PIN FIXATION BOLT SHORT
KY-0006
ROUND PIN FIXATION BOLT LONG



KY-0007
SPECIAL PIN
FIXATION BOLT 2



KY-0800
WIRE
FIXATION BUCLE



SE-7010
BUSING



HALFPIN SUPPORTING CUBES

KY-7101 1 HOLE
KY-7102 2 HOLES
KY-7103 3 HOLES
KY-7104 4 HOLES
KY-7105 5 HOLES



KY-0004
SPECIAL
FIXATION BOLT 1



KY-7001
CUBE SET SCREW



SE-0020
THREADED SOCKED L 20 MM
SE-0040
THREADED SOCKED L 40 MM




KE-0113
MULTIPLE ANGLED
UNIVERSAL HALFHINGE



KE-0112
UNIVERSAL
HALFHINGE



KY-8000
UNIVERSAL
BALL JOINT



PS-0700
M6 LONG NUT



PS-0706
SERRATED WASHER



PS-0702
SLOTTED BOLT



PS-0703
SLOTTED WASHER



PS-0701
CANNULATED BOLD



PS-0212 M6 BOLT 12 MM
PS-0213 M6 BOLT 13 MM
PS-0216 M6 BOLT 16 MM
PS-0230 M6 BOLT 30 MM



PS-0012
Ø 12 MM WASHER



PS-0017
SPECIAL WASHER



PS-0600
M6 NUT



PS-0800
M6 NUT WITH FIBRE



PS-0704
CONICAL WASHER



PS-0705
CONICAL COUPLE

EXTERNAL CIRCULAR FIXATION SYSTEM (ADULT)



OLIVER WIRES

WO-1530 Ø 1,5 MM 300 MM
 WO-1537 Ø 1,5 MM 370 MM
 WO-1540 Ø 1,5 MM 400 MM
 WO-1830 Ø 1,8 MM 300 MM
 WO-1837 Ø 1,8 MM 370 MM
 WO-1840 Ø 1,8 MM 400 MM
 WO-2030 Ø 2,0 MM 300 MM
 WO-2037 Ø 2,0 MM 370 MM
 WO-2040 Ø 2,0 MM 400 MM



BAYONET POINTED WIRES

WB-1530 Ø 1,5 MM 300 MM
 WB-1537 Ø 1,5 MM 370 MM
 WB-1540 Ø 1,5 MM 400 MM
 WB-1830 Ø 1,8 MM 300 MM
 WB-1837 Ø 1,8 MM 370 MM
 WB-1840 Ø 1,8 MM 400 MM
 WB-2030 Ø 2,0 MM 300 MM
 WB-2037 Ø 2,0 MM 370 MM
 WB-2040 Ø 2,0 MM 400 MM



SCHANZ SCREWS Ø 6

FŞ-6110S	L110 MM	FŞ-6110ST	L110 MM
FŞ-6120S	L120 MM	FŞ-6120ST	L120 MM
FŞ-6130S	L130 MM	FŞ-6130ST	L130 MM
FŞ-6140S	L140 MM	FŞ-6140ST	L140 MM
FŞ-6150S	L150 MM	FŞ-6150ST	L150 MM
FŞ-6160S	L160 MM	FŞ-6160ST	L160 MM
FŞ-6170S	L170 MM	FŞ-6170ST	L170 MM
FŞ-6180S	L180 MM	FŞ-6180ST	L180 MM
FŞ-6200S	L200 MM	FŞ-6200ST	L200 MM
FŞ-6220S	L220 MM	FŞ-6220ST	L220 MM
FŞ-6240S	L240 MM	FŞ-6240ST	L240 MM
FŞ-7110S	L110 MM	FŞ-7110ST	L110 MM
FŞ-7120S	L120 MM	FŞ-7120ST	L120 MM
FŞ-7130S	L130 MM	FŞ-7130ST	L130 MM
FŞ-7140S	L140 MM	FŞ-7140ST	L140 MM
FŞ-7150S	L150 MM	FŞ-7150ST	L150 MM
FŞ-7160S	L160 MM	FŞ-7160ST	L160 MM
FŞ-7170S	L170 MM	FŞ-7170ST	L170 MM
FŞ-7180S	L180 MM	FŞ-7180ST	L180 MM
FŞ-7200S	L200 MM	FŞ-7200ST	L200 MM
FŞ-7240S	L240 MM	FŞ-7240ST	L240 MM



TROCAR POINTED WIRES

WT-1530 Ø 1,5 MM 300 MM
 WT-1537 Ø 1,5 MM 370 MM
 WT-1540 Ø 1,5 MM 400 MM
 WT-1830 Ø 1,8 MM 300 MM
 WT-1837 Ø 1,8 MM 370 MM
 WT-1530 Ø 1,8 MM 400 MM
 WT-2030 Ø 2,0 MM 300 MM
 WT-2037 Ø 2,0 MM 370 MM
 WT-2040 Ø 2,0 MM 400 MM
 WT-0830 Ø 0,8 MM 300 MM
 WT-1030 Ø 1,0 MM 300 MM
 WT-1230 Ø 1,2 MM 300 MM
 WT-2530 Ø 2,5 MM 300 MM
 WT-3030 Ø 3,0 MM 300 MM



FULLY THREADED HALF PIN Ø 6

FŞT - 6120 L120 MM
 FŞT - 6150 L150 MM
 FŞT - 6180 L180 MM
 FŞT - 6200 L200 MM
 FŞT - 6220 L220 MM
 FŞT - 6120T L120 MM
 FŞT - 6150T L150 MM
 FŞT - 6180T L180 MM
 FŞT - 6200T L200 MM
 FŞT - 6220T L220 MM



SLOTTED SCHANZ SCREWS

WP-4150 Ø 4 MM x 150 MM
 WP-4180 Ø 4 MM x 180 MM
 WP-5150 Ø 5 MM x 150 MM
 WP-5180 Ø 5 MM x 180 MM
 WP-6150 Ø 6 MM x 150 MM
 WP-6180 Ø 6 MM x 180 MM

EXTERNAL CIRCULAR FIXATION SYSTEM (ADULT)



BONE (BONIT) COATED HALF PINS STAINLESS STEEL

BONE (BONIT) COATED HALF PINS TITANIUM

TAMF-4030S Ø 4 MM L30 MM	TAMF-6060S Ø 6 MM L60 MM
TAMF-4035S Ø 4 MM L35 MM	TAMF-6070S Ø 6 MM L70 MM
TAMF-4040S Ø 4 MM L40 MM	TAMF-6080S Ø 6 MM L80 MM
TAMF-4045S Ø 4 MM L45 MM	TAMF-6090S Ø 6 MM L90 MM
TAMF-4050S Ø 4 MM L50 MM	TAMF-6110S Ø 6 MM L110 MM
TAMF-5030S Ø 5 MM L30 MM	TAMF-6120S Ø 6 MM L120 MM
TAMF-5035S Ø 5 MM L35 MM	TAMF-6130S Ø 6 MM L130 MM
TAMF-5040S Ø 5 MM L40 MM	TAMF-6140S Ø 6 MM L140 MM
TAMF-5050S Ø 5 MM L50 MM	TAMF-6150S Ø 6 MM L150 MM
TAMF-5080S Ø 5 MM L80 MM	TAMF-6160S Ø 6 MM L160 MM
TAMF-5090S Ø 5 MM L90 MM	TAMF-6170S Ø 6 MM L170 MM
TAMF-6030S Ø 6 MM L30 MM	TAMF-6180S Ø 6 MM L180 MM
TAMF-6035S Ø 6 MM L35 MM	TAMF-6190S Ø 6 MM L190 MM
TAMF-6040S Ø 6 MM L40 MM	TAMF-6200S Ø 6 MM L200 MM
TAMF-6050S Ø 6 MM L50 MM	TAMF-6220S Ø 6 MM L220 MM
	TAMF-6250S Ø 6 MM L250 MM

TAMF-4030ST Ø 4 MM L30 MM	TAMF-6060ST Ø 6 MM L60 MM
TAMF-4035ST Ø 4 MM L35 MM	TAMF-6070ST Ø 6 MM L70 MM
TAMF-4040ST Ø 4 MM L40 MM	TAMF-6080ST Ø 6 MM L80 MM
TAMF-4045ST Ø 4 MM L45 MM	TAMF-6090ST Ø 6 MM L90 MM
TAMF-4050ST Ø 4 MM L50 MM	TAMF-6110ST Ø 6 MM L110 MM
TAMF-5030ST Ø 5 MM L30 MM	TAMF-6120ST Ø 6 MM L120 MM
TAMF-5035ST Ø 5 MM L35 MM	TAMF-6130ST Ø 6 MM L130 MM
TAMF-5040ST Ø 5 MM L40 MM	TAMF-6140ST Ø 6 MM L140 MM
TAMF-5050ST Ø 5 MM L50 MM	TAMF-6150ST Ø 6 MM L150 MM
TAMF-5080ST Ø 5 MM L80 MM	TAMF-6160ST Ø 6 MM L160 MM
TAMF-5090ST Ø 5 MM L90 MM	TAMF-6170ST Ø 6 MM L170 MM
TAMF-6030ST Ø 6 MM L30 MM	TAMF-6180ST Ø 6 MM L180 MM
TAMF-6035ST Ø 6 MM L35 MM	TAMF-6190ST Ø 6 MM L190 MM
TAMF-6040ST Ø 6 MM L40 MM	TAMF-6200ST Ø 6 MM L200 MM
TAMF-6050ST Ø 6 MM L50 MM	TAMF-6220ST Ø 6 MM L220 MM
	TAMF-6250ST Ø 6 MM L250 MM




COMPRESSION THREADED HALF PINS STAINLESS STEEL

COMPRESSION THREADED HALF PINS TITANIUM

FŞ-4140C Ø 4 MM L140 MM
FŞ-4150C Ø 4 MM L150 MM
FŞ-4160C Ø 4 MM L160 MM
FŞ-4170C Ø 4 MM L170 MM
FŞ-4180C Ø 4 MM L180 MM
FŞ-4200C Ø 4 MM L200 MM
FŞ-4220C Ø 4 MM L220 MM
FŞ-5150C Ø 5 MM L150 MM
FŞ-5160C Ø 5 MM L160 MM
FŞ-5170C Ø 5 MM L170 MM
FŞ-5180C Ø 5 MM L180 MM
FŞ-5200C Ø 5 MM L200 MM
FŞ-5220C Ø 5 MM L220 MM
FŞ-5250C Ø 5 MM L250 MM
FŞ-5260C Ø 5 MM L260 MM
FŞ-6150C Ø 6 MM L150 MM
FŞ-6160C Ø 6 MM L160 MM
FŞ-6170C Ø 6 MM L170 MM
FŞ-6180C Ø 6 MM L180 MM
FŞ-6200C Ø 6 MM L200 MM
FŞ-6220C Ø 6 MM L220 MM
FŞ-6250C Ø 6 MM L250 MM
FŞ-6260C Ø 6 MM L260 MM

FŞ-4140CT Ø 4 MM L140 MM
FŞ-4150CT Ø 4 MM L150 MM
FŞ-4160CT Ø 4 MM L160 MM
FŞ-4170CT Ø 4 MM L170 MM
FŞ-4180CT Ø 4 MM L180 MM
FŞ-4200CT Ø 4 MM L200 MM
FŞ-4220CT Ø 4 MM L220 MM
FŞ-5150CT Ø 5 MM L150 MM
FŞ-5160CT Ø 5 MM L160 MM
FŞ-5170CT Ø 5 MM L170 MM
FŞ-5180CT Ø 5 MM L180 MM
FŞ-5200CT Ø 5 MM L200 MM
FŞ-5220CT Ø 5 MM L220 MM
FŞ-5250CT Ø 5 MM L250 MM
FŞ-5260CT Ø 5 MM L260 MM
FŞ-6150CT Ø 6 MM L150 MM
FŞ-6160CT Ø 6 MM L160 MM
FŞ-6170CT Ø 6 MM L170 MM
FŞ-6180CT Ø 6 MM L180 MM
FŞ-6200CT Ø 6 MM L200 MM
FŞ-6220CT Ø 6 MM L220 MM
FŞ-6250CT Ø 6 MM L250 MM
FŞ-6260CT Ø 6 MM L260 MM

EXTERNAL CIRCULAR FIXATION SYSTEM (ADULT)



INOX SCHEDULE



CUBE SET SCREW L DRIVER



CUBE SET SCREW DRIVER




DRILL GUIDE



WIRE TENSIONER ROUND TYPE



DRILL



UNIVERSAL CHUCK WITH HANDLE



WIRE TENSIONER



OSTEOTOM



WIRE CUTTER




T HANDLE



L WRENCH




HAMMER




WRENCH WITH CIRCULAR END

AA 10 MM
AA 14 MM
AA 17 MM



LEFT RIGHT OPEN WRENCH



tasarımmed[®]

TASARIMMED TIBBİ MAMÜLLER SAN. ve TİC. LTD. ŞTİ.

Demirkapı, Topçular Cad. Set Üstü No.3 (34055) Eyüp - İSTANBUL
Tel. : +90 212 674 22 44 - 45 Fax : +90 212 674 22 46
web : www.tasarimmed.com e-mail : info@tasarimmed.com